cellsciences.com

PLAU

Native Human HMW Urokinase Plasminogen Activator, FITC

Catalog No.CRU113BQuantity:100 μg

CRU113C 1 mg

Alternate Names: U-plasminogen activator, uPA

Description: Urokinase-type Plasminogen Activator is a serine protease involved in degradation of the

extracellular matrix and possibly tumor cell migration and proliferation. A specific polymorphism in this gene may be associated with late-onset Alzheimer's disease and also with decreased affinity for fibrin-binding. This protein converts Plasminogen to Plasmin by specific cleavage of an Arg-Val bond in Plasminogen. Plasmin in turn cleaves this protein at a Lys-Ile bond to form a two-chain derivative in which a single disulfide bond connects the amino-terminal A-chain to the catalytically active, carboxy-terminal B-chain. This two-chain derivative is also called HMW-uPA (High Molecular Weight uPA). HMW-uPA can be further processed into LMW-uPA (Low Molecular Weight uPA) by cleavage of chain A into a short chain A (A1) and an ATF-uPA (Amino-Terminal

Fragment uPA). LMW-uPA is proteolytically active but does not bind to the uPA Receptor

(uPAR).

UniProt ID: P00749

Gene ID: 5328

Conjugate: Fluorescein isothiocyanate, FITC

Source: Human urine

Molecular Weight: 54.0 kDa

Formulation: 0.05 M Na Acetate, 0.1 M NaCl, 1 mM EDTA, 50% Glycerol, pH 5.0

Purity: 95% by SDS-PAGE analysis

Extinction Coefficient: $E^{0.1\%}_{280nm} = 1.36$

Concentration: ≥ 1.0 mg/ml, lot specific

Storage & Stability: Protect from light. Store at -80°C for up to 1 year. Upon initial thaw, prepare working

aliquots and store at -80°C. Avoid repeated freeze-thaw cycles.

Certification: Prepared from donors negative for HbsAG, anti-HCV, anti-HBc, anti-HIV 1 & 2 by FDA

Toll Free: 888-769-1246

Phone: 978-572-1070

Fax: 978-992-0298

E-mail: info@cellsciences.com

Website: <u>www.cellsciences.com</u>

approved tests.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.