

Recombinant HAV/Hepatitis A Virus P2C-P3A (aa 1392-1521)

Catalog No. CSI15709A Quantity: 100 μg

CSI15709B 0.5 mg CSI15709C 1.0 mg

Description: Forty-two antigenic domains were identified across the hepatitis A virus (HAV)

polyprotein by using a set of 237 overlapping 20-mer synthetic peptides spanning the entire HAV polyprotein. Nineteen antigenic domains were found within the structural proteins, and 22 were found within the nonstructural proteins, with 1 domain spanning the

junction of VP1 and P2A proteins. Five of these domains were considered immunodominant, as judged by both the breadth and the strength of their

immunoreactivity. One domain is located within the VP2 protein at position 57-90 aa. A second domain, located at position 767-842 aa, contains the C-terminal part of the VP1 protein and the entire P2A protein. A third domain, located at position 1403-1456 aa, comprises the C-terminal part of the P2C protein and the N-terminal half of the P3A protein. The fourth domain, located at position 1500-1519 aa, includes almost the entire P3B, and the last domain, located at position 1719-1764 aa, contains the C-terminal region of the P3C protein and the N-terminal region of the P3D protein. Four of the five most immunoreactive domains are derived from small HAV proteins and/or encompass protein cleavage sites separating different HAV proteins.

The E.Coli derived recombinant protein contains the P2C-P3A immunodominant regions,

amino acids 1392-1521

Source: E. coli

Formulation: 10 mM CBB, pH9.6 + 0.1% SDS +50% glycerol.

Purity: HAV P2C-P3A protein is >90% pure as determined by 10% PAGE (coomassie staining).

Purification Method: HAV P2C-P3A protein was purified by proprietary chromatographic technique.

Specific Activity: Immunoreactive with sera HAV-infected individuals

Storage & Stability: HAV P2C-P3A protein although stable at 4°C for 1 week, should be stored below -18°C.

Please prevent freeze thaw cycles.

Applications: HAV P2C-P3A antigen is suitable for ELISA and Western blots, excellent antigen for

Toll Free: 888-769-1246

Phone: 781-828-0610

Fax: 781-828-0542

E-mail: info@cellsciences.com

Website: www.cellsciences.com

detection of HAV with minimal specificity problems.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.