

Human Serum Off the Clot, Pooled Donors

Catalog No.	CSI20429A CSI20429B	Quantity:	50 ml 1 liter
Description:	Off the clot serum is produced without anticoagulants added, allowing whole blood to coagulate naturally and then separating the layers.		
Source:	Pooled human donors		
Appearance:	Yellow to orange liquid		
Preservatives:	None		
Country of Origin:	USA		
Storage & Stability:	Store at -20°C to -80°C for 1 year. Upon initial thawing, aliquot and store at -20°C to -80°C. Avoid repeated freeze/thaw cycles.		
Certification:	Non-reactive for HIV-1, HCV, HBV by NAT, HBsAg, HCV Ab, HIV 1 and 2 Ab, and RPR by currently approved FDA methods. However, because no test method can offer complete assurance that infectious agents are absent, this material should be handled at the Bio-safety level 2 (BSL 2) as recommended for any potentially infectious human serum or blood specimen in the CDC/NIH manual "Bio-safety in microbiological and Biometical Laboratories", 1999.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.