

MONOCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. 2ZHSP1B

Anti Human HSP90 β

Product type	Primary antibodies
Host	Mouse
Source	Ascites
Form	Hybridoma cell line (Hyb-K3701) was grown in ascitic fluid BALB/c Mouse Liquid Ammonium sulfate precipitation Physiological salt solution containing 0.1% Sodium Azide
Volume	0.4mL
Concentration	1mg/mL
Specificity	This antibody specifically reacts with Human HSP 90 β (Amino acid residues 185-335)
Antigen	<i>E.coli</i> recombinant human Heat-Shock Protein (HSP) 90 β antigen
Clone	Hyb-K3701 The splenic lymphocytes from BALB/c mouse immunized with <i>E.coli</i> recombinant human Heat-Shock Protein (HSP) 90 β were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-K3701) were screened.
Isotype	IgM

Application notes WB, ELISA, IHC
Recommended use

Recommended dilutions

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Storage Maintain at 2-8 °C (-20 °C for long term storage) in undiluted aliquots.
Avoid repeated freeze/thaw cycles.

References 1) Nemoto T, et al: J Biol Chem 272: 26179-26187, 1997.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

モノクローナル抗体

For research use only. Not for clinical diagnosis.

Catalog No. 2ZHSP1B

抗ヒト HSP90 β

Product type	一次抗体
Host	マウス
Source	BALB/c マウス腹水 ハイブリドーマ (Hyb-K3701) を BALB/c マウスに接種して得られた腹水
Form	液状 硫酸アンモニウム塩析精製 0.1%アジ化ナトリウムを含む生理的食塩水溶液
Volume	0.4mL
Concentration	1mg/mL
Specificity	ヒト HSP90 β (アミノ酸配列番号 ; 185-335) に特異的に反応
Antigen	大腸菌発現リコンビナントーヒト Heat-Shock Protein (HSP) 90 β 抗原
Clone	Hyb-K3701 大腸菌発現リコンビナントーヒト Heat-Shock Protein (HSP) 90 β を免疫した BALB/c マウスの脾臓細胞とマウスミエローマ細胞 (NS-1) を融合して得られたハイブリドーマ
Isotype	IgM

Application notes WB, ELISA, IHC
Recommended use

Recommended dilutions

Optimal dilutions/concentrations should be determined by the end user.

Staining Pattern

Cross reactivity

Storage 2~8°C 3 ヶ月間
-20°C 3 年間

(凍結融解の繰り返しは避けてください)

References 1) Nemoto T, et al: J Biol Chem 272: 26179-26187, 1997.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619