

Capto MMC ImpRes in ReadyToProcess columns

Product Information

Cat#No# Ca-279P

Product Overview

These ReadyToProcess chromatography columns are prepacked with Capto MMC ImpRes chromatography resin and are excellent high-resolution polishing of mAbs and other biomolecules. ReadyToProcess columns are validated high-performance bioprocessing columns that are supplied prepacked and ready for use.

Characteristic

High yields achieved through the high-resolution beads and selectivity of the ligand.
Efficient removal of aggregates, viruses, and main contaminants in MAb processes.
Enables use of a platform approach to MAb process development.
Displays a broad pH/salt operational window.
Suitable for polishing of antibody fragments.

Applications

For manufacturing of biopharmaceuticals for clinical phase I and II studies. Depending on the scale of operations, they can also be used for full-scale manufacturing, as well as for preclinical studies.

Maximum operating pressure

3 bar (0.3 MPa, 44 psi)

Storage

4°C to 30°C, 20% ethanol

Pack size

1 L

Dimensions

80 × 200 mm

Column volume

1 L

Capto MMC ImpRes in ReadyToProcess columns

Column i.d.

80 mm
