

Product Information

Endotoxin-Free Plasmid Maxi Kit

CAT#: EP-007

CAT#	PRODUCT	Package	Storage
EP-007	Endotoxin-Free Plasmid Maxi Kit	2 preparations	Multiple Temperatures
EP-007	Endotoxin-Free Plasmid Maxi Kit	10 preparations	Multiple Temperatures
EP-007	Endotoxin-Free Plasmid Maxi Kit	25 preparations	Multiple Temperatures

Product Information

Application	For fast and efficient purification of endotoxin free plasmid DNA from 150-200 mL of E. coli culture.
Kit Components	Buffer A1 Buffer B1 Buffer N3 Buffer RET DNA Wash Buffer Endofree Elution Buffer EzBind™ Columns RNase A
Storage	Multiple Temperatures
Shipping Condition	RT
Usage	For Research Use Only! Not For Use in Humans.