

HiScale SEC columns, packed on demand

Product Information

Cat#No#	Hi-316P
----------------	---------

Product Overview

If you like our new-generation Superdex Increase and Superose Increase size exclusion chromatography (SEC) columns, but need to purify larger volumes of up to 8.5 mL, we can pack the resin for you in HiScale columns.

Description

HiScale columns are packed on demand for you with Superdex 200 Increase, Superdex 75 Increase, and Superose 6 Increase SEC (gel filtration) resins. The columns are available in two sizes: 16 × 400 mm and 26 × 400 mm. These two column sizes enable purification of sample volumes up to 3.2 and 8.5 mL, respectively.

Characteristic

High-resolution SEC columns—small particle size and narrow particle size distribution of the packed resin ensures highly pure protein. HiScale columns packed with SEC resins on demand, please request a quote.

Storage

4°C to 30°C, 20% ethanol

Pack size

1

Dimensions

16 mm × 400 mm

Column volume

80 mL