

## Untagged column package

### Product Information

**Cat#No#** Un-278P

### Product Overview

The untagged column package combines cation and anion exchange with size exclusion chromatography for purification of most untagged or native proteins.

### Characteristic

Package contains HiTrap Capto Q ImpRes 1 × 1 mL and Capto SP ImpRes 1 × 1 mL ion exchange (IEX) columns for effective first-step purification.

You obtain final high purity with HiPrep 16/60 Sephacryl S-200 HR — a size exclusion chromatography (SEC) column with a broad fractionation range of Mr 5000 to ~ 250 000 for globular proteins.

Includes instruction with guidance and tips to get started quickly.

### Applications

Purification of most untagged or native proteins.

### Storage

4°C to 30°C

### Pack size

1 kit

### Kit contents

HiTrap Capto Q ImpRes (1 × 1 mL)

HiTrap Capto SP ImpRes (1 × 1 mL)

HiPrep 16/60 Sephacryl S-200 HR (1 × 120 mL)