



## Heat Aggregated Gamma Globulin (HAGG)

Cat. No.: CTR-006

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

**Product Overview** Heat Aggregated Gamma Globulin acts as a potent classical complement pathway activator, functioning similarly to naturally occurring immune complexes.

### SPECIFICATIONS

<b>Application</b>	HAGG will fully activate the classical complement pathway of the Normal Human Serum (NHS) within 30 minutes at 37°C. (10 µl HAGG +1.0 mL NHS)
<b>Composition</b>	Heat Aggregated Gamma Globulin (HAGG) in Phosphate Buffered Saline buffer (pH 7.3, 20 mM CaCl <sub>2</sub> and MgCl <sub>2</sub> ).
<b>Form</b>	Liquid
<b>Concentration</b>	50.0 mg/mL
<b>Storage</b>	Storage at -70°C