



Exosome Isolation Kit (For Other Body Fluids)

Cat. No.: EIK-12

This product is for research use only and is not intended for diagnostic use.

Size	30 reactions
Storage Conditions	RT/2-8°C
Shipment Conditions	Gel Pack
Sample Type	Other Body Fluids
Application	Exosome isolation and purification from Other Body Fluids
Product Introduction	This product is a rapid and convenient method for the efficient isolation and purification of exosomes from other body fluids, including cerebrospinal fluid, amniotic fluid, breast milk, saliva, alveolar lavage fluid and nasal lavage fluid. The kit produces highly pure exosomes using filtration method as compared to exosome precipitation methods.
Package Contents	
User Supplied Reagents and Equipment	1×PBS buffer (sterile)
Note!	<ul style="list-style-type: none"> • Cover the bottle cap immediately after each use to prevent evaporation, replace plastic tubes with glass ones; • This kit cannot isolate vesicles bigger than 300 nm; • The maximum volume of each reaction is 2 mL; • Do not exceed the recommended sample size, which may result in stratification and column blockage.