



## Cynomolgus Monkey Tonsil

Cat. No.: NHP-TI093

**This product is for research use only and is not intended for diagnostic use.**

**yes** Creative Biolabs provides this tissue block isolated from Cynomolgus monkeys' tonsil. The sample has been excluded from multiple virus infections including Herpes-B Virus, as well as SRV, SIV, and STLV-1, and can be applied in a series of *in vitro* test and assay experiments such as PCR, WB, immunoprecipitation, immunofluorescent flow cytometry.

### Detailed Product Description

|                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Source</b>       | Healthy cynomolgus monkey                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Applications</b> | Given their robust resemblances to humans across physiological, behavioral, immunological, and genetic aspects, non-human primates serve as crucial models in a broad range of biomedical research, especially in bridging translational research from small animal models to humans. Creative Biolabs offers both standard formats and personalized non-human primate tissue preparations, catering to diverse <i>in vitro</i> testing applications. |
| <b>System</b>       | Immune System                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Organ/Tissue</b> | Tonsil                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Shipping</b>     | Dry ice                                                                                                                                                                                                                                                                                                                                                                                                                                               |

### Technical Specifications

|                             |                                                                                                                          |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------|
| <b>Preservation Methods</b> | Snap frozen                                                                                                              |
| <b>Quality Assurance</b>    | Tissue blocks are prepared by histologists with years of experience to be sure of excellent morphology and high quality. |
| <b>Packaging</b>            | Securely packaged to maintain the tissue quality during shipping.                                                        |