

# DiagPoly™ Biotin Coated Fluorescent Polystyrene Particles, Pink, 5.0-5.9 μm

Cat.No: DNM-L102

#### **DESCRIPTION**

Description Biotin Coated Fluorescent PS Particles, Pink, 5.0-5.9 μm, 0.1% w/v, dispersed in DI

water, analytical standard.

### **APPLICATION**

Application Notes Cell sorting, immunoprecipitation, biosensing, in vivo and in vitro sensor, coating

materials for contrast agents and anti-cancer drugs

#### PRODUCT INFORMATION

**Diameter** 5.0-5.9 μm

**Ligand** Biotin

Concentration 0.1% w/v

### STORAGE AND SHIPPING

Storage Buffer Supplied in bottle, dispersed in DI water

**Storage** 2-8°C in the dark. Do not freeze.

**Shelf Life** When stored as specified the product is stable for at least 6 months.

## **TECHNOTES**

45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-633-6938 Fax: 1-631-938-8221

Email: info@cd-bioparticles.com

© CD Bioparticles All Rights Reserved



**TECHNOTES** 

To achieve optimum particle suspension, resuspend by vortexing before use.