

# DiagPoly™ Plain Fluorescent Polystyrene Particles, CyBlue, 10.0-14.0 µm

Cat.No: DNM-L118

## DESCRIPTION

|                    |                                                                                                             |
|--------------------|-------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Plain Fluorescent PS Particles, CyBlue, 10.0-14.0 µm, 1.0% w/v, dispersed in DI water, analytical standard. |
|--------------------|-------------------------------------------------------------------------------------------------------------|

## APPLICATION

|                          |                                                                                                       |
|--------------------------|-------------------------------------------------------------------------------------------------------|
| <b>Application Notes</b> | Suitable for uses in flow cytometry, fluorescence microscopy, phagocytosis studies and cell labeling. |
|--------------------------|-------------------------------------------------------------------------------------------------------|

## PRODUCT INFORMATION

|                 |              |
|-----------------|--------------|
| <b>Diameter</b> | 10.0-14.0 µm |
|-----------------|--------------|

|                         |       |
|-------------------------|-------|
| <b>Functional Group</b> | Plain |
|-------------------------|-------|

|                      |          |
|----------------------|----------|
| <b>Concentration</b> | 1.0% w/v |
|----------------------|----------|

## STORAGE AND SHIPPING

|                       |          |
|-----------------------|----------|
| <b>Storage Buffer</b> | DI water |
|-----------------------|----------|

|                |                                   |
|----------------|-----------------------------------|
| <b>Storage</b> | 2-8°C in the dark. Do not freeze. |
|----------------|-----------------------------------|

## TECHNOTES

|                  |                                                                            |
|------------------|----------------------------------------------------------------------------|
| <b>TECHNOTES</b> | To achieve optimum particle suspension, resuspend by vortexing before use. |
|------------------|----------------------------------------------------------------------------|