🕜 Tel: +1-301-363-4651 🖂 Email: cusabio@cusabio.com 🥑 Website: www.cusabio.com

## Recombinant Human Proactivator polypeptide(PSAP),partial


Product Code	CSB-EP018836HU
Relevance	Saposin-A and saposin-C stimulate the hydrolysis of glucosylceramide by beta- glucosylceramidase (EC 3.2.1.45) and galactosylceramide by beta- galactosylceramidase (EC 3.2.1.46). Saposin-C apparently acts by combining with the enzyme and acidic lipid to form an activated complex, rather than by solubilizing the substrate.Saposin-B stimulates the hydrolysis of galacto- cerebroside sulfate by arylsulfatase A (EC 3.1.6.8), GM1 gangliosides by beta- galactosidase (EC 3.2.1.23) and globotriaosylceramide by alpha-galactosidase A (EC 3.2.1.22). Saposin-B forms a solubilizing complex with the substrates of the sphingolipid hydrolases.Saposin-D is a specific sphingomyelin phosphodiesterase activator (EC 3.1.4.12).Prosaposin: Behaves as a myelinotrophic and neurotrophic factor, these effects are mediated by its G- protein-coupled receptors, GPR37 and GPR37L1, undergoing ligand-mediated internalization followed by ERK phosphorylation signaling.
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07602
Storage Buffer	Tris-based buffer,50% glycerol
Alias	Proactivator polypeptide
Product Type	Recombinant Protein
Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	SDVYCEVCEFLVKEVTKLIDNNKTEKEILDAFDKMCSKLPKSLSEECQEVVDTY GSSILSILLEEVSPELVCSMLHLCSGT
Research Area	Metabolism
Source	E.coli
Gene Names	PSAP
Expression Region	311-391aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	36.1kDa
Protein Description	Partial

1



🕜 Tel: +1-301-363-4651 🛛 Email: cusabio@cusabio.com 🥥 Website: www.cusabio.com 🧉

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.