





## NABP1 Antibody, Biotin conjugated

Product Code	CSB-PA850257HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AH0
Immunogen	Recombinant Human SOSS complex subunit B2 protein (1-204AA)
Raised In	Rabbit
Species Reactivity	Human
<b>Tested Applications</b>	ELISA
Relevance	Component of the SOSS complex, a multiprotein complex that functions downstream of the MRN complex to promote DNA repair and G2/M checkpoint. In the SOSS complex, acts as a sensor of single-stranded DNA that binds to single-stranded DNA, in particular to polypyrimidines. The SOSS complex associates with DNA lesions and influences diverse endpoints in the cellular DNA damage response including cell-cycle checkpoint activation, recombinational repair and maintenance of genomic stability. Required for efficient homologous recombination-dependent repair of double-strand breaks
	(DSBs) and ATM-dependent signaling pathways.
Form	(DSBs) and ATM-dependent signaling pathways.  Liquid
Form Conjugate	
	Liquid
Conjugate	Liquid Biotin Preservative: 0.03% Proclin 300
Conjugate Storage Buffer	Liquid Biotin Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Conjugate Storage Buffer Purification Method	Liquid Biotin  Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  >95%, Protein G purified
Conjugate Storage Buffer Purification Method Isotype	Liquid  Biotin  Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  >95%, Protein G purified  IgG
Conjugate Storage Buffer  Purification Method Isotype Clonality	Liquid  Biotin  Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  >95%, Protein G purified  IgG  Polyclonal  SOSS complex subunit B2 (Nucleic acid-binding protein 1) (Oligonucleotide/oligosaccharide-binding fold-containing protein 2A) (Sensor of single-strand DNA complex subunit B2) (Sensor of ssDNA subunit B2) (SOSS-
Conjugate Storage Buffer  Purification Method Isotype Clonality Alias	Liquid  Biotin  Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  >95%, Protein G purified  IgG  Polyclonal  SOSS complex subunit B2 (Nucleic acid-binding protein 1) (Oligonucleotide/oligosaccharide-binding fold-containing protein 2A) (Sensor of single-strand DNA complex subunit B2) (Sensor of ssDNA subunit B2) (SOSS-B2) (Single-stranded DNA-binding protein 2) (hSSB2), NABP1, OBFC2A SSB2