




Recombinant Murine polyomavirus Major capsid protein VP1

Product Code	CSB-EP326154MKKc7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P24595
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Murine polyomavirus (strain Kilham) (MPyV) (Murine pneumotropic virus)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	APTVKKRTSQNQGLSPQKSQNSVVGGIQVLDVRTGPDSITQIEAFLNPRMGK PVDSDFYGFSDNITVSADYTQDMPRIKELPCYSMAKISLPMLNEDMTCDTILM WEAISCKTEVVGVSLLTNCHSAVKRLYDNEGAGFPVQGLNFHFFSVGGGEALD LQWLWKNYRCNYPAGVAALQAAPKAAQVLDPKLKAKLTADGKFPIEAWSPDP AKNENTRYFGTYTGGLQTPPVLLQITNTTTTILLNENGVGPLCKGDGLYLASADI VGFRTQQNNKMHLRGLPRYFSIHLRKGCANPYPVSSLLNTFSSEMMPLNSWM LQVEEVRIYDGVRLPGDPDMIRYRIIWPGRLLSLIFPAMRHKHLYFFVMQAFIV L
Lead Time	3-7 business days
Research Area	others
Source	E.coli
Gene Names	N/A
Expression Region	2-373aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	42.4 kDa
Protein Description	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.