




Recombinant Enterobacteria phage PRD1 Spike protein P5 (V)

Product Code	CSB-EP326327EDO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P22536
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage PRD1 (Bacteriophage PRD1)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	ANQQIGGSTVTYNGAIPMGGPVAINSVIEIAGTEVLVDLKL DYATGKISGVQTL YI DLRDFLGDVTVTMPDTGQRITARAGTQGYYPVLSTNLMKFIVSATIDGKFPMN FINFPIALGVWPSGIKGDKGDPGAPGPAGGTVVVEDSGASFGESLLDTTSEPG KILVKRISGGSGITVTDYGDQVEIEASGGGGGGGGVTDALSLMYSTSTGGPASI AANALTD FDL SGALT VNSVGTGLTKSAAGIQLAAGKSGLYQITMTVKNNTVTTG NYLLRVKYGSSDFV VACPASSLTAGGTISLLIYCNVLGVPSLDVLKFSLCNDGA ALSNIINITA AKIN
Lead Time	Delivery time may differ from different purchasing way or location, please kindly consult your local distributors for specific delivery time.
Research Area	Others
Source	E.coli
Gene Names	V
Expression Region	2-340aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	40.3 kDa
Protein Description	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20^o/-80^o. Our default final concentration of glycerol is 50%. Customers could use it as reference.