# NKMAXBIO We support you, we believe in your research

# Recombinant human MYL6B protein

Catalog Number: ATGP1206

#### PRODUCT INFORMATION

## **Expression system**

E.coli

#### **Domain**

1-208aa

#### **UniProt No.**

P14649

#### **NCBI Accession No.**

NP 002466

#### **Alternative Names**

Myosin light chain 6B, MLC1SA

## **PRODUCT SPECIFICATION**

## **Molecular Weight**

25.2 kDa (231aa) confirmed by MALDI-TOF

#### Concentration

1mg/ml (determined by Bradford assay)

#### **Formulation**

Liquid in. 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 2mM DTT

#### **Purity**

> 95% by SDS-PAGE

#### Tag

His-Tag

## **Application**

SDS-PAGE

# **Storage Condition**

Can be stored at +2C to +8C for 1 week. For long term storage, aliquot and store at -20C to -80C. Avoid repeated freezing and thawing cycles.

#### **BACKGROUND**

#### **Description**

MYL6B is a heavy chain regulator found in smooth muscle and non-muscle Myosin complexes. Contractile activity in smooth muscle is regulated by the calcium/calmodulin-dependent phosphorylation of Myosin light chain by Myosin light chain kinase. MYL6B does not bind calcium during contraction. It is primarily found in a hexamer consisting of four light chains and two heavy chains. It most commonly interacts with Myosin Va, an Actin based motor that can move in large steps. MYL6B is expressed in most tissues with neurons and smooth muscle tissue having the highest expression. Recombinant human MYL6B protein, fused to His-tag at N-



# NKMAXBio We support you, we believe in your research

# Recombinant human MYL6B protein

Catalog Number: ATGP1206

terminus, was expressed in E. coli and purified by using conventional chromatography techniques.

# **Amino acid Sequence**


MGSSHHHHHH SSGLVPRGSH MGSMPPKKDV PVKKPAGPSI SKPAAKPAAA GAPPAKTKAE PAVPQAPQKT QEPPVDLSKV VIEFNKDQLE EFKEAFELFD RVGDGKILYS QCGDVMRALG QNPTNAEVLK VLGNPKSDEL KSRRVDFETF LPMLQAVAKN RGQGTYEDYL EGFRVFDKEG NGKVMGAELR HVLTTLGEKM TEEEVETVLA GHEDSNGCIN YEAFLKHILS V

#### **General References**

Mitsui T., et al. (1992) J Biol Chem. 267(23):16727-35. Takeya K., et al. (2008) Am J Physiol Renal Physiol. 294(6):1487-92.

# **DATA**

#### **SDS-PAGE**



3ug by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

