

## TNF-alpha / TNFA Antibody (APC), Rabbit MAb

Catalog Number: 10602-R101-A



EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human TNF-alpha protein (Catalog#10602-HNAE)
Reagents:	APC-conjugated Rabbit monoclonal antibody
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human TNF-alpha / TNFA (rh TNF-alpha / TNFA; Catalog#10602-HNAE; NP_000585.2; Val77-Leu233) and conjugated with APC under optimum conditions, the unreacted APC was removed.
Ig Type:	Rabbit IgG
Clone ID:	101
Specificity:	Human TNF-alpha / TNFA
Concentration:	5 μl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.1% sodium azide
Storage:	This antibody can be stored at $2^{\circ}\text{C-8}^{\circ}\text{C}$ for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Alternative Names:	TNF
APPLICATIONS	
Applications:	FCM
RECOMMENDED CONCENTRATION	

Please Note: Optimal concentrations/dilutions should be determined by the end user.




## TNF-alpha / TNFA Antibody (APC), Rabbit MAb

Catalog Number: 10602-R101-A



EliteRmab® is a registered trademark of Sino Biological Inc.



Flow cytometric analysis of Human TNF-α expression on human peripheral blood lymphocytes. Human peripheral blood mononulcear cells were stimulated for 4-6 hours with PMA and lonomycin in the presence of GolgiPlug. The cells were treated according to manufacturer's manual (BD Pharmingen™ Cat. No. 554714), stained with APC-conjugated anti-Human TNF-α and FITC-conjugated anti-Human CD3. The dot plots were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.