

# CD200 Antibody, Mouse MAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 10886-MM17

## GENERAL INFORMATION

<b>Immunogen:</b>	Recombinant Human CD200 Protein (Catalog#10886-H08H)
<b>Preparation</b>	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD200 (rh CD200; Catalog#10886-H08H; NP_005935.4; Met1-Gly232). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
<b>Ig Type:</b>	Mouse IgG1
<b>Clone ID:</b>	17
<b>Specificity:</b>	Human CD200
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## APPLICATIONS

**Applications:** ELISA,FCM

## RECOMMENDED CONCENTRATION

<b>Flow Cytometry</b>	FCM: 1:25-1:100
<b>ELISA</b>	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human CD200 / OX2.


**Please Note:** *Optimal concentrations/dilutions should be determined by the end user.*

# CD200 Antibody, Mouse MAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 10886-MM17



Flow cytometric analysis of Human CD200 expression on human whole blood lymphocytes. Cells were stained with purified anti-Human CD200, then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to [www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html](http://www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html) for technical protocols.