

# FYB Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological  
Biological Solution Specialist

Catalog Number: 203556-T44

## GENERAL INFORMATION

Immunogen:	E. coli-derived Human FYB fragment
Preparation	Produced in rabbits immunized with E. coli-derived Human FYB fragment, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human FYB
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

## APPLICATIONS

Applications:	WB, IHC-P, IP
---------------	---------------

## RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:50-1:200
Western Blot	WB: 1:500-1:2000
Immunoprecipitation	IP:1-5μL/mg of lysate


**Please Note: Optimal concentrations/dilutions should be determined by the end user.**

# FYB Antibody, Rabbit PAb, Antigen Affinity Purified


Catalog Number: 203556-T44



Sino Biological  
Biological Solution Specialist



Immunohistochemical staining of human FYB in human lymph node with rabbit polyclonal antibody at 1:100 dilution, formalin-fixed paraffin embedded sections.



Anti-FYB rabbit polyclonal antibody at 1:500 dilution

Lane A: TF-1 Whole Cell Lysate

Lane B: Jurkat Whole Cell Lysate


Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.  
Performed under reducing conditions.

Predicted band size: 85 kDa  
Observed band size: 120 kDa  
(We are unsure as to the identity of these extra bands.)



FYB was immunoprecipitated using:

Lane A: 0.5 mg Jurkat Whole Cell Lysate

4 µL anti-FYB rabbit polyclonal antibody and 60 µg of Immunomagnetic beads Protein A/G.

Primary antibody:

Anti-FYB rabbit polyclonal antibody, at 1:100 dilution

Secondary antibody:

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution

Developed using the ECL technique.  
Performed under reducing conditions.

Predicted band size: 85 kDa  
Observed band size: 120 kDa