# Rhesus DC-SIGN / CD209 Protein(His Tag)

Catalog Number: 90319-C07H



## **General Information**

#### Gene Name Synonym:

CD209: DC-SIGN

#### **Protein Construction:**

A DNA sequence encoding the rhesus CD209 (AAK74185.1) (Lys62-Glu381) was expressed with a polyhistidine tag at the N-terminus.

Source: Rhesus

Expression Host: HEK293 Cells

**QC** Testing

**Purity:** > 95 % as determined by SDS-PAGE.

**Endotoxin:** 

< 1.0 EU per µg protein as determined by the LAL method.

Stability:

Samples are stable for up to twelve months from date of receipt at -70  $^{\circ}$ C

Predicted N terminal: His

**Molecular Mass:** 

The recombinant rhesus CD209 consists of 339 amino acids and predicts a molecular mass of 38.8 kDa.

### Formulation:

Lyophilized from sterile PBS, pH 7.4.

Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Specific concentrations are included in the hardcopy of COA. Please contact us for any concerns or special requirements.

# **Usage Guide**

### Storage:


Store it under sterile conditions at  $\text{-}20\,^\circ\!\!\!\mathrm{C}$  to  $\text{-}80\,^\circ\!\!\!\mathrm{C}$  upon receiving. Recommend to aliquot the protein into smaller quantities for optimal storage.

Avoid repeated freeze-thaw cycles.

#### Reconstitution:

Detailed reconstitution instructions are sent along with the products.

#### SDS-PAGE:



# **Protein Description**

Dendritic cell (DC)-specific intercellular adhesion molecule 3 (ICAM-3) grabbing nonintegrin (DC-SIGN), also known as CD209, is a type II transmembrane protein on DCs with a C-type lectin extracellular domain, is capable of binding ICAM-3 on resting T cells in the secondary lymphoid organs, providing the initial contact between these cells during the establishment of cell-mediated immunity. It is not only a pattern recognition receptor but implicated in immunoregulation of DCs. It has important role in mediating DC adhesion, migration, inflammation, activating primary T cell, triggering immune response and participating in immune escape of pathogens and tumors. DC-SIGN also mediates capture and internalization of viral, bacterial, and fungal pathogens by dendritic cells, such as HIV-1, Ebola virus, cytomegalovirus, Dengue virus, and hepatitis C virus. DC-SIGN is unique in that it regulates adhesion processes, such as DC trafficking and T-cell synapse formation, as well as antigen capture. Moreover, even though several C-type lectins have been shown to bind HIV-1, DC-SIGN does not only capture HIV-1 but also protects it in early endosomes allowing HIV-1 transport by DC to lymphoid tissues, where it enhances trans infection of T cells.

#### References

1.Geijtenbeek TB, et al. (2002) DC-SIGN, a C-type lectin on dendritic cells that unveils many aspects of dendritic cell biology. J Leukoc Biol. 71(6): 921-31. 2.Masso M. (2003) DC-SIGN points the way to a novel mechanism for HIV-1 transmission. MedGenMed. 5(2): 2. 3.Zhou T, et al. (2006) DC-SIGN and immunoregulation. Cell Mol Immunol. 3(4): 279-83.

Manufactured By Sino Biological Inc., FOR RESEARCH USE ONLY. NOT FOR USE IN HUMANS.

For US Customer: Fax: 267-657-0217 • Tel: 215-583-7898

Global Customer: Fax :+86-10-5862-8288 
■ Tel:+86-400-890-9989 
■ http://www.sinobiological.com