# Rhesus Interferon Gamma/IFN gamma/IFNG Gene ORF cDNA clone in cloning vector

Catalog Number: CG90008-M



#### **General Information**

Gene: interferon-gamma

Official Symbol: IFNG

Synonym: IFNG

Source: Rhesus

cDNA Size: 498bp

**RefSeq:** NM\_001032905.1

Plasmid: pMD-monIFNG

#### **Description**

Lot: Please refer to the label on the tube

#### **Sequence Description:**

Identical with the Gene Bank Ref. ID sequence. Please check the sequence information before order.

Vector:

pMD18-T Simple

#### **Quality control:**

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

#### Sequencing primer list:

M13-47: 5' GCCAGGGTTTTCCCAGTCACGAC 3'

RV-M: 5' GAGCGGATAACAATTTCACACAGG3'

Other M13 primers can also be used as sequencing primers.

#### Shipping carrier:

Each tube contains approximately 10 µg of lyophilized plasmid.

#### Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

#### **Plasmid Resuspension protocol**

- 1. Centrifuge at 5,000 × g for 5 min.
- 2. Carefully open the tube and add 100  $\mu l$  of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than  $5000\times g$ .
- 5. Store the plasmid at  $-20 \,^{\circ}$ C.

#### The plasmid is ready for:

- · Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

## *E.coli* strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 $\alpha$  and TOP10F $^{\prime}$ .

For U.S. Customer: Fax: 267-657-0217 Tel: 215-583-7898

Non-U.S. Customer: Fax: +86-10-5862-8288 Tel: +86-400-890-9989 <a href="http://www.sinobiological.com">http://www.sinobiological.com</a>

### Rhesus Interferon Gamma/IFN gamma/IFNG Gene ORF cDNA clone in cloning vector

**Catalog Number:** CG90008-M




#### **Vector Information**

pMD18-T Simple Vector is a high-efficiency TA cloning vector constructed from pUC18, of which the initial multiple cloning sites (MCS) were destroyed. The pMD18-T Simple Vector is 2.6kb in size and contains the ampicillin resistance gene for selection. The coding sequence was inserted by TA cloning at site 425.

Notes: The direction of cDNA insertion into the TA-cloning vector is random, maybe forward or reverse. For insert orientation information, please feel free to contact us.

#### Physical Map of pMD18-T Simple (MCS destroyed):



Sequencing Primer M13-47

5'-CGCCAGGGTTTTCCCAGTCACGACGTTGTAAAACGACGGCCAGTGCC

EcoR V AAAGAAGCATGACGGCAAGTGGAC**GATATC**TCCAGAGGATCGCCGGGAA Cloning Site:425 ...gtggacgatT atctccaga... (TA cloning Site) ...cacctgcta Ttagaggtct...

CCGAGGACGAGTTCGTAATCATGGTCATAGCTGTTT<u>CCTGTGTGAAATTGTT</u> Sequencing Primer RV-M ATCCGCTC -3'

Non-U.S. Customer: Fax: +86-10-5862-8288 Tel: +86-400-890-9989 http://www.sinobiological.com