Rhesus CD40/TNFRSF5 Gene ORF cDNA clone in cloning vector

Catalog Number: CG90097-G



General Information

Gene: CD40 molecule, TNF receptor

superfamily member 5

Official Symbol: CD40

Synonym: CD40

Source: Rhesus

cDNA Size: 837bp

RefSeq: XM_001104333.2

Plasmid: pGEM-cynoCD40

Description

Lot: Please refer to the label on the tube

Sequence Description:

Identical with the Gene Bank Ref. ID sequence. Please check the sequence information before order.

Vector:

pGEM-T

Quality control:

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list:

M13-47: 5' GCCAGGGTTTTCCCAGTCACGAC 3'

RV-M: 5' GAGCGGATAACAATTTCACACAGG3'

Other M13 primers can also be used as sequencing primers.

Shipping carrier:

Each tube contains approximately 10 µg of lyophilized plasmid.

Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

Plasmid Resuspension protocol

- 1. Centrifuge at 5,000 × g for 5 min.
- 2. Carefully open the tube and add 100 μl of sterile water to dissolve the DNA.
- Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 5000 \times g.
- 5. Store the plasmid at $-20 \,^{\circ}$ C.

The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 α and TOP10F $^{\prime}$.

For U.S. Customer: Fax: 267-657-0217 Tel: 215-583-7898

Non-U.S. Customer: Fax: +86-10-5862-8288 Tel: +86-400-890-9989 http://www.sinobiological.com

Rhesus CD40/TNFRSF5 Gene ORF cDNA clone in cloning vector

Catalog Number: CG90097-G




Vector Information

The pGEM-T vector is a high-efficiency TA cloning vector which contains multiple cloning sites as shown below. The pGEM-T vector is 3.0kb in size and contains the ampicillin resistance gene for selection. The coding sequence was inserted by TA cloning.

Notes: The direction of cDNA insertion into the TA-cloning vector is random, maybe forward or reverse. For insert orientation information, please feel free to contact us.

Physical Map of pGEM-T:



For U.S. Customer: Fax: 267-657-0217 Tel: 215-583-7898

Non-U.S. Customer: Fax: +86-10-5862-8288 Tel: +86-400-890-9989 http://www.sinobiological.com