# Mouse IL-23 (IL23A & IL12B Heterodimer) Protein

Catalog Number: CT028-M08H



# **General Information**

#### Gene Name Synonym:

IL-23; p19

#### **Protein Construction:**

A DNA sequence encoding the mouse IL23A (Q9EQ14) (Met1-Ala196) was fused with a polyhistidine tag at the C-terminus, constructed the plasmid 1; A DNA sequence encoding mouse IL12B (P43432) (Met1-Ser335) was fused with a polyhistidine tag at the C-terminus, constructed the plasmid 2. The two plasmids were co-expressed and the mouse IL23 heterodimer was purified.

Source: Mouse

Expression Host: HEK293 Cells

**QC** Testing

Purity: > 95 % as determined by SDS-PAGE

# **Bio Activity:**

1. Immobilized Human IL23 Receptor hFc(Cat:13840-H02H) at 2  $\mu$ g/ml (100  $\mu$ l/well) can bind Mouse IL-23 His(Cat:CT028-M08H), the EC<sub>50</sub> of Mouse IL-23 His is 50-250 ng/mL.

2. Measured by its ability to induce IL17 secretion by mouse splenocytes. The  $\rm ED_{50}$  for this effect is 0.5-4 ng/mL.

#### **Endotoxin:**

< 1.0 EU per µg of the protein as determined by the LAL method

Predicted N terminal: Val 22 & Met 23

#### **Molecular Mass:**

The recombinant mouse IL23 (mouse IL23A&IL12B Heterodimer) comprises 510 (186+324) amino acids and has a calculated molecular mass of 58.3 (21.1+37.2) KDa. The apparent molecular mass of the recombinant protein is approximately 47,45 and 26 KDa respectively in SDS-PAGE under reducing conditions.

#### Formulation:

Lyophilized from sterile PBS, pH 7.4.

Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Specific concentrations are included in the hardcopy of COA. Please contact us for any concerns or special requirements.

# **Usage Guide**

## Stability & Storage:

Samples are stable for twelve months from date of receipt at -20°C to -80°C.

Store it under sterile conditions at -20°C to -80°C upon receiving. Recommend to aliquot the protein into smaller quantities for optimal storage.

# Avoid repeated freeze-thaw cycles.

# Reconstitution:

Detailed reconstitution instructions are sent along with the products.

#### SDS-PAGE:

