Human E-Cadherin ORF mammalian expression plasmid, C-Flag tag



Catalog Number: HG10204-CF

General Information

Gene: cadherin 1, type 1, E-cadherin (epithelial)

(CDH1)

Official Symbol: CDH1

Synonym: CDH1, UVO, CDHE, ECAD, LCAM, Arc-

1, CD324

Source: Human

cDNA Size: 2649bp

RefSeq: NM_004360.3

Plasmid: pCMV3-CDH1-Flag

Description

Lot: Please refer to the label on the tube

Sequence Description:

Identical with the Gene Bank Ref. ID sequence except for the point mutation 2076 T/C not causing the amino acid variation.

Restriction site: HindIII + XbaI (6kb + 2.7kb)

Vector: pCMV3-C-FLAG

Shipping carrier:

Each tube contains approximately 10 µg of lyophilized plasmid.

Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

Quality control :

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list:

pCMV3-F: 5' CAGGTGTCCACTCCCAGGTCCAAG 3'
pcDNA3-R: 5' GGCAACTAGAAGGCACAGTCGAGG 3'

Or

Forward T7: 5' TAATACGACTCACTATAGGG 3'

ReverseBGH: 5' TAGAAGGCACAGTCGAGG 3'

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

Plasmid Resuspension protocol

- 1. Centrifuge at $5,000 \times g$ for 5 min.
- 2.Carefully open the tube and add 100 μl of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4.Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than $5000\times g$.
- 5.Store the plasmid at -20 °C.

The plasmid is ready for:

- · Restriction enzyme digestion
- · PCR amplification
- · E. coli transformation
- DNA sequencing

E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 α and TOP10F $^{\prime}$.

Human E-Cadherin ORF mammalian expression plasmid, C-Flag tag



Catalog Number: HG10204-CF

Vector Information

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

- •Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.
- Hygromycin resistance gene for selection of mammalian cell lines.
- A Kozak consensus sequence to enhance mammalian expression.

Vector Name pCMV3-C-FLAG

Vector Size 6158bp

Vector Type Mammalian Expression Vector Expression Method Constitutive, Stable / Transient

Promoter CMV
Antibiotic Resistance Kanamycin

Selection In Mammalian Cells Hygromycin

Protein Tag FLAG

Physical Map of Plasmid:

