

# Human HDAC4 Gene ORF cDNA clone expression plasmid, C-His tag



Sino Biological  
Biological Solution Specialist

Catalog Number: HG11807-CH

## General Information

**Gene :** histone deacetylase 4  
**Official Symbol :** HDAC4  
**Synonym :** AHO3; BDMR; HA6116; HD4; HDAC-4; HDAC-A; HDACA  
**Source :** Human  
**cDNA Size:** 3300bp (cDNA Size= Gene + linker +Tags)  
**RefSeq :** NM\_006037.3  
**Plasmid:** pCMV3-HDAC4-His

## Description

**Lot :** Please refer to the label on the tube

### Sequence Description :

Identical with the Gene Bank Ref. ID sequence except for the point mutations: 1269G/A, 2565C/T not causing the amino acid variation.

**Restriction site:** KpnI(two restriction sites) + XbaI(6kb+2.87kb+0.44kb)

**Vector :** pCMV3-C-His

### Quality control :

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

### Sequencing primer list :

|            |                                |
|------------|--------------------------------|
| pCMV3-F:   | 5' CAGGTGTCCACTCCCAGGTCCAAG 3' |
| pcDNA3-R : | 5' GGCAACTAGAAGGCACAGTCGAGG 3' |
| Or         |                                |

|              |                            |
|--------------|----------------------------|
| Forward T7 : | 5' TAATACGACTCACTATAGGG 3' |
| ReverseBGH : | 5' TAGAAGGCACAGTCGAGG 3'   |

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

### Shipping carrier :

Each tube contains approximately 10 µg of lyophilized plasmid.

### Storage :

The lyophilized plasmid can be stored at ambient temperature for three months.

## Plasmid Resuspension protocol

1. Centrifuge at 5,000×g for 5 min.
2. Carefully open the tube and add 100 µl of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 5000×g.
5. Store the plasmid at -20 °C.

### The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- *E. coli* transformation
- DNA sequencing

## *E.coli* strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5α and TOP10F'.

Manufactured By Sino Biological Inc., FOR RESEARCH USE ONLY. NOT FOR USE IN HUMANS.

For U.S. Customer: Fax: 267-657-0217 Tel: 215-583-7898

Non-U.S. Customer: Fax: +86-10-5862-8288 Tel: +86-400-890-9989

<http://www.sinobiological.com>

# Human HDAC4 Gene ORF cDNA clone expression plasmid, C-His tag



Sino Biological  
Biological Solution Specialist

Catalog Number: HG11807-CH


## Vector Information

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

- Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.
- Hygromycin resistance gene for selection of mammalian cell lines.
- A Kozak consensus sequence to enhance mammalian expression.

|                      |                                  |
|----------------------|----------------------------------|
| Vector name          | pCMV3-C-His                      |
| Vector size          | 6164bp                           |
| Vector Type          | Mammalian Expression Vector      |
| Expression Method    | Constitutive ,Stable / Transient |
| Promoter             | CMV                              |
| Bacterial Resistance | Kanamycin                        |
| Selection In Cells   | Hygromycin                       |
| Protein tag          | His                              |

## Physical Map of Plasmid :



Manufactured By Sino Biological Inc., FOR RESEARCH USE ONLY. NOT FOR USE IN HUMANS.

For U.S. Customer: Fax: 267-657-0217 Tel: 215-583-7898

Non-U.S. Customer: Fax: +86-10-5862-8288 Tel: +86-400-890-9989

<http://www.sinobiological.com>