# Human GALNT15 Gene cDNA clone plasmid

Catalog Number: HG13692-G



### **General Information**

Gene: UDP-N-acetyl-alpha-D-

> galactosamine:polypeptide N-

acetylgalactosaminyltransferase 15

Official Symbol: **GALNT15** 

Synonym: PIH5, GALNT7, GALNT13, GALNT15,

GALNTL2

Source: Human

cDNA Size: 1920bp

RefSeq: BC014789

Plasmid: pGEM-GALNTL2

**Description** 

Lot: Please refer to the label on the tube

**Sequence Description:** 

Identical with the Gene Bank Ref. ID sequence.

Vector:

pGEM-T

Shipping carrier:

Each tube contains approximately 10 µg of lyophilized plasmid.

Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

**Quality control:** 

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list:

M13-47: 5' GCCAGGGTTTTCCCAGTCACGAC 3'

RV-M: 5' GAGCGGATAACAATTTCACACAGG 3'

Other M13 primers can also be used as sequencing primers.

### **Plasmid Resuspension protocol**

- 1. Centrifuge at  $5,000 \times g$  for 5 min.
- 2.Carefully open the tube and add 100 µl of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 5000×g.
- 5.Store the plasmid at -20 °C.

### The plasmid is ready for:

- · Restriction enzyme digestion
- · PCR amplification
- · E. coli transformation
- DNA sequencing

## E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5α and TOP10F'.

# **Human GALNT15 Gene cDNA clone plasmid**

Catalog Number: HG13692-G



### **Vector Information**

The pGEM-T vector is a high-efficiency TA cloning vector which contains multiple cloning sites as shown below. The pGEM-T vector is 3.0kb in size and contains the amplicin resistance gene for selection. The coding sequence was inserted by TA cloning.

### Physical Map of pGEM-T:

