# Human DNASE1 Gene cDNA clone plasmid

Catalog Number: HG13801-G



#### **General Information**

Gene: deoxyribonuclease I

Official Symbol: **DNASE1** 

Synonym: DKFZp686H0155, DRNI, DNL1,

FLJ38093, FLJ44902, DNASE1

Source: Human

cDNA Size: 849bp

RefSeq: BC029437

Plasmid: pGEM-DNASE1

### **Description**

Lot: Please refer to the label on the tube

#### **Sequence Description:**

Identical with the Gene Bank Ref. ID sequence.

Vector:

pGEM-T

#### Shipping carrier:

Each tube contains approximately 10 µg of lyophilized plasmid.

### Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

#### **Quality control:**

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

#### Sequencing primer list:

M13-47: 5' GCCAGGGTTTTCCCAGTCACGAC 3'

RV-M: 5' GAGCGGATAACAATTTCACACAGG 3'

Other M13 primers can also be used as sequencing primers.

### **Plasmid Resuspension protocol**

- 1. Centrifuge at 5,000 × g for 5 min.
- 2.Carefully open the tube and add 100 µl of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 5000×g.
- 5.Store the plasmid at -20 °C.

### The plasmid is ready for:

- · Restriction enzyme digestion
- · PCR amplification
- · E. coli transformation
- DNA sequencing

## E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5α and TOP10F'.

# **Human DNASE1 Gene cDNA clone plasmid**

Sino Biological Inc.

**Biological Solution Specialist** 

**Catalog Number:** HG13801-G

#### **Vector Information**

The pGEM-T vector is a high-efficiency TA cloning vector which contains multiple cloning sites as shown below. The pGEM-T vector is 3.0kb in size and contains the amplicin resistance gene for selection. The coding sequence was inserted by TA cloning.

#### Physical Map of pGEM-T:

