Human COMMD10 Gene cDNA clone plasmid

Catalog Number: HG14008-G

General Information

Gene : COMM domain containing 10

Official Symbol :	COMMD10
onnoiai oynnoor .	

- PTD002, COMMD10 Synonym :
- Source : Human
- cDNA Size: 609bp
- BC005179 **RefSeq:**
- Plasmid: pGEM-COMMD10

Description

Lot : Please refer to the label on the tube

Sequence Description :

Identical with the Gene Bank Ref. ID sequence.

Vector :

pGEM-T

Shipping carrier :

Each tube contains approximately 10 µg of lyophilized plasmid.

Storage :

The lyophilized plasmid can be stored at ambient temperature for three months.

Quality control:

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list :

M13-47 :	5' GCCAGGGTTTTCCCAGTCACGAC 3'

RV-M: 5' GAGCGGATAACAATTTCACACAGG 3'

Other M13 primers can also be used as sequencing primers.



Plasmid Resuspension protocol

1.Centrifuge at $5,000 \times g$ for 5 min.

2.Carefully open the tube and add 100 μ l of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature. 4.Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than $5000 \times g$. 5.Store the plasmid at -20 °C.

The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 α and TOP10F'.

Human COMMD10 Gene cDNA clone plasmid


Catalog Number: HG14008-G



Vector Information

The pGEM-T vector is a high-efficiency TA cloning vector which contains multiple cloning sites as shown below. The pGEM-T vector is 3.0kb in size and contains the amplicin resistance gene for selection. The coding sequence was inserted by TA cloning.

Physical Map of pGEM-T :



Manufactured By Sino Biological Inc., FOR RESEARCH USE ONLY. NOT FOR USE IN HUMANS. Fax :+86-10-51029969 • Tel:+86- 400-890-9989 • <u>http://www.sinobiological.com</u>