# Human PFK1/PFKM Gene ORF cDNA clone expression plasmid, N-HA tag



**Catalog Number:** HG14133-NY

**General Information** 

Gene: phosphofructokinase, muscle

Official Symbol: **PFKM** 

ATP-PFK; GSD7; PFK-1; PFK1; PFKA; Synonym:

PFKX; PPP1R122

Source: Human

cDNA Size: 2385bp

RefSeq: BC021203

Plasmid: pCMV3-HA-PFKM

Description

Please refer to the label on the tube Lot:

**Sequence Description:** 

Identical with the Gene Bank Ref. ID sequence.

Restriction site: KpnI + XbaI(6kb+2.39kb)

Vector: pCMV3-N-HA

Quality control:

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list:

pCMV3-F: 5' CAGGTGTCCACTCCCAGGTCCAAG 3'

pcDNA3-R: 5' GGCAACTAGAAGGCACAGTCGAGG 3'

Or

Forward T7: 5' TAATACGACTCACTATAGGG 3'

ReverseBGH: 5' TAGAAGGCACAGTCGAGG 3'

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

#### Shipping carrier:

Each tube contains approximately 10  $\mu g$  of lyophilized plasmid.

#### Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

## **Plasmid Resuspension protocol**

- 1. Centrifuge at  $5,000 \times g$  for 5 min.
- 2.Carefully open the tube and add 100  $\mu l$  of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than  $5000 \times g$ .
- 5.Store the plasmid at -20 °C.

#### The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- · E. coli transformation
- · DNA sequencing

# E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5α and TOP10F'.

# Human PFK1/PFKM Gene ORF cDNA clone expression plasmid, N-HA tag



Catalog Number: HG14133-NY

#### **Vector Information**

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

- •Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.
- Hygromycin resistance gene for selection of mammalian cell lines.
- A Kozak consensus sequence to enhance mammalian expression.

Vector Name pCMV3-N-HA

Vector Size 6101bp

Vector Type Mammalian Expression Vector Expression Method Constitutive, Stable / Transient

Promoter CMV

Antibiotic Resistance Kanamycin Selection In Mammalian Cells Hygromycin

Protein Tag HA

### **Physical Map of Plasmid:**

