Human VAMP2 Gene ORF cDNA clone expression plasmid, C-His tag



Catalog Number: HG14320-CH

		Shi
General Information		
Gene :	vesicle-associated membrane protein 2 (synaptobrevin 2)	Eac Sto
Official Symbol :	VAMP2	The
Synonym :	SYB2; VAMP-2	thre
Source :	Human	Pla
cDNA Size:	396bp (cDNA Size= Gene + linker +Tags)	
RefSeq :	BC002737	1. C
Plasmid:	pCMV3-VAMP2-His	2. C
Description		
Lot : Please refer to the label on the tube		3. C
		tem
Sequence Description :		4. E
Identical with the Gene Bank Ref. ID sequence.		th a
Restriction site:	KpnI + Xbal(6kb+0.4kb)	the
Vector :	pCMV3-C-His	5. S
Quality control :		
The plasmid is confirmed by full-length sequencing with primers		

in the sequencing primer list. Sequencing primer list :

pCMV3-F:	5' CAGGTGTCCACTCCCAGGTCCAAG 3'
pcDNA3-R :	5' GGCAACTAGAAGGCACAGTCGAGG 3'
Or	
Forward T7 :	5' TAATACGACTCACTATAGGG 3'
ReverseBGH :	5' TAGAAGGCACAGTCGAGG 3'

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

Shipping	carrier :
----------	-----------

Each tube contains approximately 10 µg of lyophilized plasmid. Storage :

The lyophilized plasmid can be stored at ambient temperature for three months.

Plasmid Resuspension protocol

1. Centrifuge at $5,000 \times g$ for 5 min.

2. Carefully open the tube and add 100 μ l of sterile water to

dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room

temperature.

4. Briefly vortex the tube and then do a quick spin to concentrate

the liquid at the bottom. Speed is less than $5000 \times g$.

5. Store the plasmid at -20 °C.

The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 α and TOP10F'.

Human VAMP2 Gene ORF cDNA clone expression plasmid, C-His tag



Catalog Number: HG14320-CH

Vector Information

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

•Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.

• Hygromycin resistance gene for selection of mammalian cell lines.

• A Kozak consensus sequence to enhance mammalian expression.

Vector name	pCMV3-C-His
Vector size	6164bp
Vector Type	Mammalian Expression Vector
Expression Method	Constiutive ,Stable / Transient
Promoter	CMV
Bacterial Resistance	Kanamycin
Selection In Cells	Hygromycin
Protein tag	His

Physical Map of Plasmid :

