# Homo sapiens PMP22 cDNA Clone

Catalog Number: HG14519-G-F



#### **General Information**

Gene: peripheral myelin protein 22

Official Symbol: PMP22

Synonym: DSS, HNPP, CMT1A, CMT1E, GAS-3,

Sp110, HMSNIA

Source: Homo sapiens

cDNA Size: 483

RefSeq: BC019040

**Description** 

Lot: Please refer to the label on the tube

**Sequence Description:** 

Identical with the Gene Bank Ref. ID sequence.

Restriction site: Kpnl + Xhol

Vector:

pCMV / hygro - FLAG

Shipping carrier:

Each tube contains approximately 10 µg of lyophilized plasmid.

Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

**Quality control:** 

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list:

5' CTAGAGAACCCACTGCTTACTGGC 3' pcDNA3-L:

pcDNA3-R: 5' GGCAACTAGAAGGCACAGTCGAGG 3'

Or

Forward T7: 5' TAATACGACTCACTATAGGG 3'

ReverseBGH: 5' TAGAAGGCACAGTCGAGG 3'

pcDNA3-L and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

### **Plasmid Resuspension protocol**

- 1. Centrifuge at  $5,000 \times g$  for 5 min.
- 2. Carefully open the tube and add 100 µl of sterile water to dissolve the DNA.
- 3.Close the tube and incubate for 10 minutes at room temperature.
- 4.Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than  $5000\times q$ .
- 5. Store the plasmid at  $-20 \,^{\circ}$ C.

#### The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

## E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5α and TOP10F'.

# Homo sapiens PMP22 cDNA Clone

**Catalog Number:** HG14519-G-F




#### **Vector Information**

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

- · Human cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.
- Hygromycin resistance gene for selection of mammalian cell
- A Kozak consensus sequence to enhance mammalian expression.

Vector Name	pCMV / hygro - FLAG
Vector Size	5681bp
Vector Type	Mammalian Expression Vector
Expression Method	Constitutive, Stable / Transient
Promoter	CMV
Antibiotic Resistance	Ampicillin
Selection In Mammalian Cells	Hygromycin
Protein Tag	FLAG

### **Physical Map of Plasmid:**



http://www.sinobiological.com/Mammalian-Expression-Vectors-a-4666.html for the vector sequence.

Fax:+86-10-51029969