# **Human HEYL Gene cDNA clone plasmid**

Catalog Number: HG14720-G



### **General Information**

Gene: hairy/enhancer-of-split related with

YRPW motif-like

Official Symbol: HEYL

Synonym: HEY3, HRT3, bHLHb33

Source: Human

cDNA Size: 987bp

RefSeg: BC006087

Plasmid: pGEM-HEYL

**Description** 

Lot: Please refer to the label on the tube

**Sequence Description:** 

Identical with the Gene Bank Ref. ID sequence.

Vector:

pGEM-T

Shipping carrier:

Each tube contains approximately 10 µg of lyophilized plasmid.

Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

**Quality control:** 

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list:

M13-47: 5' GCCAGGGTTTTCCCAGTCACGAC 3'

RV-M: 5' GAGCGGATAACAATTTCACACAGG 3'

Other M13 primers can also be used as sequencing primers.

## **Plasmid Resuspension protocol**

- 1. Centrifuge at  $5,000 \times g$  for 5 min.
- 2.Carefully open the tube and add 100  $\,\mu l$  of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4.Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than  $5000 \times g$ .
- 5.Store the plasmid at -20 °C.

## The plasmid is ready for:

- · Restriction enzyme digestion
- · PCR amplification
- · E. coli transformation
- DNA sequencing

# *E.coli* strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 $\alpha$  and TOP10F $^{\prime}$ .

Fax:+86-10-51029969

● Tel:+86- 400-890-9989 ● <a href="http://www.sinobiological.com">http://www.sinobiological.com</a>

# **Human HEYL Gene cDNA clone plasmid**

Catalog Number: HG14720-G



### **Vector Information**

The pGEM-T vector is a high-efficiency TA cloning vector which contains multiple cloning sites as shown below. The pGEM-T vector is 3.0kb in size and contains the amplicin resistance gene for selection. The coding sequence was inserted by TA cloning.

### Physical Map of pGEM-T:

