# Homo sapiens CHRFAM7A cDNA Clone

Catalog Number: HG15733-G-F

# **General Information**

Gene :	CHRNA7 (cholinergic receptor, nicotinic, alpha 7, exons 5-10) and FAM7A (family with sequence similarity 7A, exons A-E) fusion
Official Symbol :	CHRFAM7A
Synonym :	D-10, CHRNA7, CHRNA7-DR1
Source :	Homo sapiens
cDNA Size:	1239
RefSeq :	BC101348

## Description

Lot : Please refer to the label on the tube

#### **Sequence Description :**

Identical with the Gene Bank Ref. ID sequence.

Restriction site: HindIII + Xhol

Vector :

pCMV / hygro - FLAG

#### Shipping carrier :

Each tube contains approximately 10 µg of lyophilized plasmid.

#### Storage :

The lyophilized plasmid can be stored at ambient temperature for three months.

#### **Quality control :**

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

#### Sequencing primer list :

pcDNA3-L :	5' CTAGAGAACCCACTGCTTACTGGC 3'
pcDNA3-R :	5' GGCAACTAGAAGGCACAGTCGAGG 3'
Or	
Forward T7 :	5' TAATACGACTCACTATAGGG 3'
ReverseBGH :	5' TAGAAGGCACAGTCGAGG 3'

pcDNA3-L and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.



## **Plasmid Resuspension protocol**

1.Centrifuge at  $5,000 \times g$  for 5 min.

2.Carefully open the tube and add 100  $\mu l$  of sterile water to dissolve the DNA.

3.Close the tube and incubate for 10 minutes at room temperature.

4.Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than  $5000 \times g$ .

5.Store the plasmid at -20  $^{\circ}$ C.

#### The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

# *E.coli* strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5  $\alpha$  and TOP10F  $\dot{}$  .

# Homo sapiens CHRFAM7A cDNA Clone

Catalog Number: HG15733-G-F

### **Vector Information**


All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

• Human cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.

• Hygromycin resistance gene for selection of mammalian cell lines.

• A Kozak consensus sequence to enhance mammalian expression.

#### **Physical Map of Plasmid :**



Please refer to http://www.sinobiological.com/Mammalian-Expression-Vectors-a-4666.html for the vector sequence.

Vector Name	pCMV / hygro - FLAG
Vector Size	5681bp
Vector Type	Mammalian Expression Vector
Expression Method	Constitutive, Stable / Transient
Promoter	CMV
Antibiotic Resistance	Ampicillin
Selection In Mammalian Cells	Hygromycin
Protein Tag	FLAG

