Human VDAC1 Gene ORF cDNA clone expression plasmid, C-Flag tag



Catalog Number: HG17755-CF

General Information		Shipping carrier :	
Gene :	voltage-dependent anion channel 1	Each tube contains approximately 10 µg of lyophilized plasmid. Storage :	
Official Symbol :	VDAC1	The lyophilized plasmid can be stored at ambient temperature for	
Synonym :	PORIN; VDAC-1	three months.	
Source :	Human		
cDNA Size:	891bp (cDNA Size= Gene + linker +Tags)	Plasmid Resuspension protocol	
RefSeq :	BC090042.1	1. Centrifuge at 5,000 $ imes$ g for 5 min.	
Plasmid:	pCMV3-VDAC1-Flag	2. Carefully open the tube and add 100 μl of sterile water to	
Description		dissolve the DNA.	
Lot : Please refer to the label on the tube		3. Close the tube and incubate for 10 minutes at room	
Sequence Description :		temperature.	
Identical with the Gene Bank Ref. ID sequence.		4. Briefly vortex the tube and then do a quick spin to concentrate	
Restriction site:	Kpnl + Xbal(6kb+0.89kb)	the liquid at the bottom. Speed is less than 5000 $ imes$ g.	
Vector :	pCMV3-C-FLAG	5. Store the plasmid at -20 °C.	
Quality control :			
The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.		The plasmid is ready for:	
Sequencing primer list :		Restriction enzyme digestion	

PCR amplification

• DNA sequencing

but not limited)

• E. coli transformation

Sequencing primer list :

pCMV3-F:	5' CAGGTGTCCACTCCCAGGTCCAAG 3'
pcDNA3-R :	5' GGCAACTAGAAGGCACAGTCGAGG 3'
Or	
Forward T7 :	5' TAATACGACTCACTATAGGG 3'
ReverseBGH :	5' TAGAAGGCACAGTCGAGG 3'

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 α and TOP10F'.

E.coli strains for transformation (recommended

Human VDAC1 Gene ORF cDNA clone expression plasmid, C-Flag tag



Catalog Number: HG17755-CF

Vector Information

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

•Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.

• Hygromycin resistance gene for selection of mammalian cell lines.

• A Kozak consensus sequence to enhance mammalian expression.

Vector name	pCMV3-C-FLAG
Vector size	6158bp
Vector Type	Mammalian Expression Vector
Expression Method	Constiutive ,Stable / Transient
Promoter	CMV
Bacterial Resistance	Kanamycin
Selection In Cells	Hygromycin
Protein tag	FLAG

Physical Map of Plasmid :

