# Human GART ORF mammalian expression plasmid, N-GFPSpark tag



Catalog Number: HG19133-ANG

**General Information** 

**Gene:** phosphoribosylglycinamide

formyltransferase,

phosphoribosylglycinamide synthetase,

phosphoribosylaminoimidazole

synthetase

Official Symbol: GART

**Synonym:** AIRS; GARS; PAIS; PGFT; PRGS;

**GARTF** 

Source: Human

cDNA Size: 1302bp

**RefSeq:** NM\_175085.2

**Description** 

**Lot:** Please refer to the label on the tube

**Vector:** pCMV3-N-GFPSpark

Shipping carrier:

Each tube contains approximately 10  $\mu g$  of lyophilized plasmid.

Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

**Quality control:** 

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list:

pCMV3-F: 5' CAGGTGTCCACTCCCAGGTCCAAG 3'
pcDNA3-R: 5' GGCAACTAGAAGGCACAGTCGAGG 3'

Or

Forward T7: 5' TAATACGACTCACTATAGGG 3'
ReverseBGH: 5' TAGAAGGCACAGTCGAGG 3'

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

### **Plasmid Resuspension protocol**

- 1. Centrifuge at 5,000×g for 5 min.
- 2. Carefully open the tube and add 100  $\mu$ l of sterile water to dissolve the DNA.
- Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than  $5000 \times g$ .
- 5. Store the plasmid at -20 °C.

#### The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

## *E.coli* strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 $\alpha$  and TOP10F  $\dot{}$  .

## Human GART ORF mammalian expression plasmid, N-GFPSpark tag



Catalog Number: HG19133-ANG

### **Vector Information**

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

•Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.

 Hygromycin resistance gene for selection of mammalian cell lines.

 A Kozak consensus sequence to enhance mammalian expression. Vector Name pCMV3-N-GFPSpark

Vector Size 6788bp

Vector Type Mammalian Expression Vector

Expression Method Constitutive, Stable / Transient

Promoter CMV

Antibiotic Resistance Kanamycin

Selection In Mammalian Cells

Hygromycin

Protein Tag

**GFPSpark**