

Product Introduction

The **Magnetic Separator-1.5 (2 tubes)** can hold two microcentrifuge tubes (1.5 mL) for immunoprecipitation (IP) applications (standard IP, co-IP, chromatin IP (ChIP), or RNA IP (RIP)). It contains high-energy magnets that the magnetic beads are pulled down in seconds.

Simplicity—the separator is used for Sino Biological immunomagnetic beads conjugated Protein A, Protein G, Protein L or IP antibodies decreasing the protocol time, making IP simple and fast.

Manufactured By Sino Biological Inc., FOR RESEARCH USE ONLY. NOT FOR USE IN HUMANS.

Tel: 215-583-7898 (US); +86-400-890-9989(Global)

<http://www.sinobiological.com>