



Product Introduction

Sinobiological 16 well Magnetic Separator was designed for IP, Co-IP and Ch-IP. It contains 16 wells (1.5ml or 2mL).

The 16 well magnetic separator was made with eight pieces of N35 magnet and PMMA plates. If you have to prepare several samples for your immunoprecipitation testes, it will be helpful. With which you can put 16 tubes a time, and the samples are enough for two pieces of gels for western blotting detection. When you prepare same sample with the magnetic separator, it can treat up to 32 mL a time.

Manufactured By Sino Biological Inc., FOR RESEARCH USE ONLY. NOT FOR USE IN HUMANS.

Tel: 215-583-7898 (US); +86-400-890-9989(Global)

<http://www.sinobiological.com>