# Mouse DNAJC30 Gene ORF cDNA clone expression plasmid, C-Flag tag



Catalog Number: MG53163-CF

**General Information** 

Gene: DnaJ (Hsp40) homolog, subfamily C,

member 30

Official Symbol: DNAJC30

**Synonym:** 1300007M11Rik; Wbscr18

Source: Mouse

cDNA Size: 699bp

**RefSeq:** NM\_025362.3

Plasmid: pCMV3-mDNAJC30-Flag

Description

Lot: Please refer to the label on the tube

**Sequence Description:** 

Identical with the Gene Bank Ref. ID sequence.

**Restriction site:** Kpnl + Xbal(6kb+0.7kb)

Vector: pCMV3-C-FLAG

**Quality control:** 

The plasmid is confirmed by full-length sequencing with primers

in the sequencing primer list.

Sequencing primer list :

pCMV3-F: 5' CAGGTGTCCACTCCCAGGTCCAAG 3'

pcDNA3-R: 5' GGCAACTAGAAGGCACAGTCGAGG 3'

Or

Forward T7: 5' TAATACGACTCACTATAGGG 3'

ReverseBGH: 5' TAGAAGGCACAGTCGAGG 3'

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

#### **Shipping carrier:**

Each tube contains approximately 10 µg of lyophilized plasmid.

### Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

## **Plasmid Resuspension protocol**

- 1. Centrifuge at 5,000×g for 5 min.
- 2. Carefully open the tube and add 100  $\mu l$  of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than  $5000 \times g$ .
- 5. Store the plasmid at  $-20 \,^{\circ}$ C.

#### The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

# *E.coli* strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 $\alpha$  and TOP10F'.

# Mouse DNAJC30 Gene ORF cDNA clone expression plasmid, C-Flag tag



Catalog Number: MG53163-CF

#### **Vector Information**

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

- •Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.
- Hygromycin resistance gene for selection of mammalian cell lines.
- A Kozak consensus sequence to enhance mammalian expression.

Vector Name pCMV3-C-FLAG

Vector Size 6158bp

Vector Type Mammalian Expression Vector Expression Method Constitutive, Stable / Transient

Promoter CMV
Antibiotic Resistance Kanamycin
Selection In Mammalian Cells Hygromycin

Protein Tag FLAG

### Physical Map of Plasmid:

