

Catalog Number: RG80880-CH

## **General Information**

| Gene : decapping enzyme, s | scavenger |
|----------------------------|-----------|
|----------------------------|-----------|

| Official Symbol : | DCPS       |
|-------------------|------------|
| Synonym :         | Hint-5     |
| Source :          | Rat        |
| cDNA Size:        | 1008bp     |
| RefSeq :          | BC105900.1 |

# Description

| Lot : | Please refer to the label on the tube |
|-------|---------------------------------------|
| LUL.  |                                       |

Vector : pCMV3-C-His

#### Shipping carrier :

Each tube contains approximately 10  $\mu g$  of lyophilized plasmid.

#### Storage :

The lyophilized plasmid can be stored at ambient temperature for three months.

#### **Quality control :**

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

#### Sequencing primer list :

| pCMV3-F:     | 5' CAGGTGTCCACTCCCAGGTCCAAG 3' |
|--------------|--------------------------------|
| pcDNA3-R :   | 5' GGCAACTAGAAGGCACAGTCGAGG 3' |
| Or           |                                |
| Forward T7 : | 5' TAATACGACTCACTATAGGG 3'     |
| ReverseBGH : | 5' TAGAAGGCACAGTCGAGG 3'       |
|              |                                |

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

# **Plasmid Resuspension protocol**

- 1. Centrifuge at  $5,000 \times g$  for 5 min.
- 2. Carefully open the tube and add 100  $\mu$ l of sterile water to dissolve the DNA.
- Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than  $5000 \times g$ .
- 5. Store the plasmid at -20  $^\circ\!\mathrm{C}.$

## The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

# *E.coli* strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5  $\alpha$  and TOP10F  $\dot{}$  .



**Catalog Number:** RG80880-CH

## **Vector Information**

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

•Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.

· Hygromycin resistance gene for selection of mammalian cell lines.

• A Kozak consensus sequence to enhance mammalian expression.

| Vector Name                  | pCMV3-C-His                      |
|------------------------------|----------------------------------|
| Vector Size                  | 6164bp                           |
| Vector Type                  | Mammalian Expression Vector      |
| Expression Method            | Constitutive, Stable / Transient |
| Promoter                     | CMV                              |
| Antibiotic Resistance        | Kanamycin                        |
| Selection In Mammalian Cells | Hygromycin                       |
| Protein Tag                  | His                              |