Rat PSMC3 ORF mammalian expression plasmid, C-Flag tag



Catalog Number: RG81004-CF

General Information

| Gene : | proteasome (prosome, subunit, ATPase 3 | macropain) | 26S |
|-------------------|--|------------|-----|
| Official Symbol : | PSMC3 | | |
| Synonym : | Psmc3 | | |

Source :

| cDNA Size: | 1329bp |
|------------|--------|
| | |

Rat

RefSeq : BC062019.1

Description

Lot : Please refer to the label on the tube

Vector : pCMV3-C-FLAG

Shipping carrier :

Each tube contains approximately 10 µg of lyophilized plasmid.

Storage :

The lyophilized plasmid can be stored at ambient temperature for three months.

Quality control :

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list :

| pCMV3-F: | 5' CAGGTGTCCACTCCCAGGTCCAAG 3' |
|--------------|--------------------------------|
| pcDNA3-R : | 5' GGCAACTAGAAGGCACAGTCGAGG 3' |
| Or | |
| Forward T7 : | 5' TAATACGACTCACTATAGGG 3' |
| ReverseBGH : | 5' TAGAAGGCACAGTCGAGG 3' |

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

Plasmid Resuspension protocol

- 1. Centrifuge at $5,000 \times g$ for 5 min.
- 2. Carefully open the tube and add 100 μ l of sterile water to dissolve the DNA.
- Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than $5000 \times g$.
- 5. Store the plasmid at -20 $^{\circ}$ C.

The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 α and TOP10F $\dot{}$.

Rat PSMC3 ORF mammalian expression plasmid, C-Flag tag



Catalog Number: RG81004-CF

Vector Information

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

•Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.

• Hygromycin resistance gene for selection of mammalian cell lines.

• A Kozak consensus sequence to enhance mammalian expression.

| Vector Name | pCMV3-C-FLAG |
|------------------------------|----------------------------------|
| Vector Size | 6158bp |
| Vector Type | Mammalian Expression Vector |
| Expression Method | Constitutive, Stable / Transient |
| Promoter | CMV |
| Antibiotic Resistance | Kanamycin |
| Selection In Mammalian Cells | Hygromycin |
| Protein Tag | FLAG |