Influenza A H1N1 (A/Brevig Mission/1/1918) Matrix protein 1 / M1 ORF mammalian expression plasmid (Codon Optimized)



Catalog Number: VG40211-UT

General Information

Gene: H1N1 (A/Brevig Mission/1/1918) Matrix

protein 1 / M1.

Official Symbol:

Synonym: Matrix protein 1, M1

Source: H₁N₁

cDNA Size: 759bp

Plasmid: pCMV3-H1N1-B-1918-M1

Description

Lot: Please refer to the label on the tube

Sequence Description:

The translated amino acid sequence is identical with Q8BAC3.1.

Restriction site: KpnI + XbaI (6.1kb + 0.76kb)

Vector: pCMV3-untagged

Shipping carrier:

Each tube contains approximately 10 µg of lyophilized plasmid.

Storage:

The lyophilized plasmid can be stored at ambient temperature

for three months.

Quality control:

The plasmid is confirmed by full-length sequencing with primers in the sequencing primer list.

Sequencing primer list:

pCMV3-F: 5' CAGGTGTCCACTCCCAGGTCCAAG 3'

pcDNA3-R: 5' GGCAACTAGAAGGCACAGTCGAGG 3'

Or

5' TAATACGACTCACTATAGGG 3' Forward T7:

ReverseBGH: 5' TAGAAGGCACAGTCGAGG 3'

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

Plasmid Resuspension protocol

- Centrifuge at $5,000 \times g$ for 5 min.
- Carefully open the tube and add 100 µl of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than 5000×g.
- 5. Store the plasmid at -20 °C.

The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- E. coli transformation
- DNA sequencing

E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5 α and TOP10F'.

Influenza A H1N1 (A/Brevig Mission/1/1918) Matrix protein 1 / M1 ORF mammalian expression plasmid (Codon Optimized)



Catalog Number: VG40211-UT

Vector Information

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

- •Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.
- Hygromycin resistance gene for selection of mammalian cell lines.
- A Kozak consensus sequence to enhance mammalian expression.

Vector Name pCMV3-untagged

Vector Size 6223bp

Vector Type Mammalian Expression Vector Expression Method Constitutive, Stable / Transient

None

Promoter CMV
Antibiotic Resistance Ampicillin
Selection In Mammalian Cells Hygromycin

Protein Tag

Physical Map of Plasmid:

