# SARS coronavirus (SARS-CoV) (strain WH20) Plpro ORF mammalian expression plasmid, N-Flag tag



**Catalog Number:** VG40524-NF

**General Information** 

Official Symbol: **PLpro** 

Synonym: **PLpro** 

Source: SARS

cDNA Size: 1002bp (cDNA Size= Gene + linker

+Tags)

RefSeq: AY772062.1(nt4895-5851)

Plasmid: pCMV3-Flag-PLpro

Description

Lot: Please refer to the label on the tube

**Sequence Description:** 

Identical with the Gene Bank Ref. ID sequence.

Restriction site: KpnI + XbaI(6kb+1kb)

Vector: pCMV3-N-FLAG

Quality control:

The plasmid is confirmed by full-length sequencing with primers

in the sequencing primer list.

Sequencing primer list:

5' CAGGTGTCCACTCCCAGGTCCAAG 3' pCMV3-F: pcDNA3-R: 5' GGCAACTAGAAGGCACAGTCGAGG 3' Or

Forward T7: 5' TAATACGACTCACTATAGGG 3'

ReverseBGH: 5' TAGAAGGCACAGTCGAGG 3'

pCMV3-F and pcDNA3-R are designed by Sino Biological Inc. Customers can order the primer pair from any oligonucleotide supplier.

#### Shipping carrier:

Each tube contains approximately 10  $\mu g$  of lyophilized plasmid.

#### Storage:

The lyophilized plasmid can be stored at ambient temperature for three months.

## **Plasmid Resuspension protocol**

- 1. Centrifuge at 5,000×g for 5 min.
- 2. Carefully open the tube and add 100  $\mu$ l of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin to concentrate the liquid at the bottom. Speed is less than  $5000 \times g$ .
- 5. Store the plasmid at  $-20 \,^{\circ}$ C.

### The plasmid is ready for:

- Restriction enzyme digestion
- PCR amplification
- · E. coli transformation
- · DNA sequencing

# E.coli strains for transformation (recommended but not limited)

Most commercially available competent cells are appropriate for the plasmid, e.g. TOP10, DH5α and TOP10F'.

Website: <a href="http://www.sinobiological.com">http://www.sinobiological.com</a>

# SARS coronavirus (SARS-CoV) (strain WH20) Plpro ORF mammalian expression plasmid, N-Flag tag



Catalog Number: VG40524-NF

## **Vector Information**

All of the pCMV vectors are designed for high-level stable and transient expression in mammalian hosts. High-level stable and non-replicative transient expression can be carried out in most mammalian cells. The vectors contain the following elements:

- •Human enhanced cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells.
- Hygromycin resistance gene for selection of mammalian cell lines.
- A Kozak consensus sequence to enhance mammalian expression.

Vector name pCMV3-N-FLAG

Vector size 6098bp

Vector Type Mammalian Expression Vector

Expression Method Constitutive ,Stable / Transient

Promoter CMV

Bacterial Resistance Kanamycin
Selection In Cells Hygromycin

Protein tag FLAG

#### **Physical Map of Plasmid:**

