ITCH / AIP4 Antibody, Mouse MAb



Catalog Number: 11131-MM01


GENERAL INFORMATION	
Immunogen:	Recombinant Human ITCH / AIP4 protein (Catalog#11131-HNCE)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human ITCH / AIP4 (rh ITCH / AIP4; Catalog#11131-HNCE; NP_113671.3; Arg526-Glu903). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	01
Specificity:	Human ITCH / AIP4
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at $2^{\circ}\text{C}-8^{\circ}\text{C}$ for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C . Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	ADMFD,AIF4,AIP4,dJ468O1.1,NAPP1
APPLICATIONS	
Applications:	FCM,ICC/IF
RECOMMENDED CONCENTRATION	
ICC/IF	ICC/IF: 1:20-1:100
Flow Cytometry	FCM: 1:25-1:100

Please Note: Optimal concentrations/dilutions should be determined by the end user.


ITCH / AIP4 Antibody, Mouse MAb

Catalog Number: 11131-MM01





Immunofluorescence staining of Human ITCH in Hela cells. Cells were fixed with 4% PFA, permeabilzed with 1% Triton X-100 in PBS, blocked with 10% serum, and incubated with Mouse anti-Human ITCH monoclonal antibody (1:60) at 37 °C 1 hour. Then cells were stained with the Alexa Fluor® 488-conjugated Goat Anti-mouse IgG secondary antibody (green) and counterstained with DAPI (blue). Positive staining was localized to cytoplasm.



Flow cytometric analysis of Human ITCH expression on HeLa cells. The cells were treated according to manufacturer's manual (BD Pharmingen™ Cat. No. 554714), stained with purified anti-Human ITCH, then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.