

Anti-IFIT2 antibody



Description Unconjugated Rabbit polyclonal to IFIT2

Model STJ190316

Host Rabbit

Reactivity Human

Applications ELISA, WB

Immunogen Synthesized peptide derived from human IFIT2 protein.

Immunogen Region 200-280aa

Gene ID <u>3433</u>

Gene Symbol <u>IFIT2</u>

Dilution range WB 1:500-2000 ELISA 1:5000-20000

Specificity IFIT2 Polyclonal Antibody detects endogenous levels of protein.

Purification IFIT2 antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Note For Research Use Only (RUO).

Protein Name Interferon-induced protein with tetratricopeptide repeats 2 IFIT-2 ISG-54 K

Interferon-induced 54 kDa protein IFI-54K P54

Molecular Weight 51 kDa

Clonality Polyclonal

Conjugation Unconjugated

Isotype IgG

Formulation Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.

Concentration 1 mg/ml

Storage Instruction Store at -20°C, and avoid repeat freeze-thaw cycles.

Database Links <u>HGNC:5409OMIM:147040</u>

Alternative Names Interferon-induced protein with tetratricopeptide repeats 2 IFIT-2 ISG-54 K

Interferon-induced 54 kDa protein IFI-54K P54

Function IFN-induced antiviral protein which inhibits expression of viral messenger

RNAs lacking 2'-O-methylation of the 5' cap. The ribose 2'-O-methylation would provide a molecular signature to distinguish between self and non-self mRNAs by the host during viral infection. Viruses evolved several ways to evade this restriction system such as encoding their own 2'-O-methylase for their mRNAs or by stealing host cap containing the 2'-O-methylation (cap snatching mechanism). Binds AU-rich viral RNAs, with or without 5' triphosphorylation, RNA-binding is required for antiviral activity. Can

promote apoptosis.

Sequence and Domain Family The C-terminal part folds into a super-helical structure and has an extensively

positively-charged nucleotide-binding channel on its inner surface.

Cellular Localization Cytoplasm Endoplasmic reticulum

St John's Laboratory Ltd

F +44 (0)207 681 2580

T +44 (0)208 223 3081 **E** info

E info@stjohnslabs.com

W http://www.stjohnslabs.com/