

## **Anti-DHRS9** antibody



**Description** Unconjugated Rabbit polyclonal to DHRS9

Model STJ190743

**Host** Rabbit

**Reactivity** Human, Mouse, Rat

**Applications** ELISA, WB

**Gene ID** 10170

Gene Symbol DHRS9

**Dilution range** WB 1:500-2000 ELISA 1:5000-20000

**Specificity** DHRS9 Polyclonal Antibody detects endogenous levels of protein.

**Tissue Specificity** Highly expressed in trachea and epidermis. Detected at lower levels in spinal

cord, bone marrow, brain, tongue, esophagus, heart, colon, testis, placenta,

lung, skeletal muscle and lymph node.

**Purification** DHRS9 antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

**Note** For Research Use Only (RUO).

Protein Name Dehydrogenase/reductase SDR family member 9 3-alpha hydroxysteroid

dehydrogenase 3-alpha-HSD NADP-dependent retinol

dehydrogenase/reductase RDH-E2 RDHL Retinol dehydrogenase 15 Short

chain dehydrogenase/reductase

Molecular Weight 35 kDa

**Clonality** Polyclonal

**Conjugation** Unconjugated

**Isotype** IgG

**Formulation** Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.

**Concentration** 1 mg/ml

**Storage Instruction** Store at -20°C, and avoid repeat freeze-thaw cycles.

Database Links <u>HGNC:16888OMIM:612131</u>

Alternative Names Dehydrogenase/reductase SDR family member 9 3-alpha hydroxysteroid

dehydrogenase 3-alpha-HSD NADP-dependent retinol

dehydrogenase/reductase RDH-E2 RDHL Retinol dehydrogenase 15 Short

chain dehydrogenase/reductase

**Function** 3-alpha-hydroxysteroid dehydrogenase that converts 3-alpha-

tetrahydroprogesterone (allopregnanolone) to dihydroxyprogesterone and 3-alpha-androstanediol to dihydroxyprogesterone. May play a role in the biosynthesis of retinoic acid from retinaldehyde, but seems to have low

activity with retinoids. Can utilize both NADH and NADPH.

Cellular Localization Microsome membrane Endoplasmic reticulum membrane. Associated with

microsomal membranes.

St John's Laboratory Ltd

**F** +44 (0)207 681 2580 **T** +44 (0)208 223 3081

W http://www.stjohnslabs.com/ E info@stjohnslabs.com