


Anti-ACSVL6 antibody



Description Rabbit polyclonal to ACSVL6.

Model STJ91460

Host Rabbit

Reactivity Human

ELISA, IF, IHC, WB **Applications**

Synthesized peptide derived from human ACSVL6 **Immunogen**

450-530 aa, Internal **Immunogen Region**

Gene ID 10998

Gene Symbol SLC27A5

WB 1:500-1:2000IHC 1:100-1:300IF 1:200-1:1000ELISA 1:40000 **Dilution range**

ACSVL6 Polyclonal Antibody detects endogenous levels of ACSVL6 protein. **Specificity**

Predominantly expressed in liver. **Tissue Specificity**

Purification The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Note For Research Use Only (RUO).

Protein Name Bile acyl-CoA synthetase BACS Bile acid-CoA ligase BA-CoA ligase BAL

Cholate--CoA ligase Fatty acid transport protein 5 FATP-5 Fatty-acid-

coenzyme A ligase, very long-chain 3 Solute carrier family 27 m

Molecular Weight 75 kDa

Polyclonal **Clonality**

Conjugation Unconjugated

Isotype IgG

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Concentration 1 mg/ml

Storage Instruction Store at -20°C, and avoid repeat freeze-thaw cycles.

Database Links <u>HGNC:10999OMIM:603314</u>

Alternative Names Bile acyl-CoA synthetase BACS Bile acid-CoA ligase BA-CoA ligase BAL

Cholate--CoA ligase Fatty acid transport protein 5 FATP-5 Fatty-acid-coenzyme A ligase, very long-chain 3 Solute carrier family 27 m

Function Acyl-CoA synthetase involved in bile acid metabolism. Proposed to catalyze

the first step in the conjugation of C24 bile acids (choloneates) to glycine and taurine before excretion into bile canaliculi by activating them to their CoA thioesters. Seems to activate secondary bile acids entering the liver from the enterohepatic circulation. In vitro, also activates 3-alpha,7-alpha,12-alpha-trihydroxy-5-beta-cholestanate (THCA), the C27 precursor of cholic acid

deriving from the de novo synthesis from cholesterol.

Cellular Localization Endoplasmic reticulum membrane

St John's Laboratory Ltd F +44 (0)207 681 2580

T +44 (0)208 223 3081

W http://www.stjohnslabs.com/ E info@stjohnslabs.com