

Anti-AKAP 220 antibody



Description Rabbit polyclonal to AKAP 220.

Model STJ91530

Host Rabbit

Reactivity Human

Applications ELISA, IF, IHC

ImmunogenSynthesized peptide derived from human AKAP 220

Immunogen Region 1730-1810 aa, C-terminal

Gene ID <u>11215</u>

Gene Symbol <u>AKAP11</u>

Dilution range IHC 1:100-1:300IF 1:200-1:1000ELISA 1:10000

Specificity AKAP 220 Polyclonal Antibody detects endogenous levels of AKAP 220

protein.

Tissue Specificity Expressed in heart, brain, lung, liver, kidney, testis and ovary. Weakly

expressed in skeletal muscle, pancreas and spleen.

Purification The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Note For Research Use Only (RUO).

Protein Name A-kinase anchor protein 11 AKAP-11 A-kinase anchor protein 220 kDa

AKAP 220 hAKAP220 Protein kinase A-anchoring protein 11 PRKA11

Molecular Weight 210.512 kDa

Clonality Polyclonal

Conjugation Unconjugated

Isotype IgG

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Concentration 1 mg/ml

Storage Instruction Store at -20°C, and avoid repeat freeze-thaw cycles.

Database Links HGNC:369OMIM:604696

Alternative Names A-kinase anchor protein 11 AKAP-11 A-kinase anchor protein 220 kDa

AKAP 220 hAKAP220 Protein kinase A-anchoring protein 11 PRKA11

Function Binds to type II regulatory subunits of protein kinase A and anchors/targets

them.

Sequence and Domain Family RII-alpha binding site, predicted to form an amphipathic helix, could

participate in protein-protein interactions with a complementary surface on the

R-subunit dimer.

Cellular Localization Cytoplasm, cytoskeleton, microtubule organizing center,

centrosome. Cytoplasmic in premeiotic pachytene spermatocytes and in the

centrosome of developing postmeiotic germ cells, while a

midpiece/centrosome localization was found in elongating spermatocytes and

mature sperm.

St John's Laboratory Ltd

F +44 (0)207 681 2580

T +44 (0)208 223 3081

W http://www.stjohnslabs.com/ E info@stjohnslabs.com