
Anti-ATP5H antibody





Description Rabbit polyclonal to ATP5H.

Model STJ91771

Host Rabbit

Reactivity Human, Mouse, Rat

Applications ELISA, IHC, WB

Immunogen Synthesized peptide derived from human ATP5H

Immunogen Region 80-160 aa, C-terminal

Gene ID <u>10476</u>

Gene Symbol ATP5H

Dilution range WB 1:500-1:2000IHC 1:100-1:300ELISA 1:20000

Specificity ATP5H Polyclonal Antibody detects endogenous levels of ATP5H protein.

Purification The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Note For Research Use Only (RUO).

Protein Name ATP synthase subunit d, mitochondrial ATPase subunit d

Molecular Weight 26 kDa

Clonality Polyclonal

Conjugation Unconjugated

Isotype IgG

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Concentration 1 mg/ml

Storage Instruction Store at -20°C, and avoid repeat freeze-thaw cycles.

Database Links <u>HGNC:845</u>OMIM:NA

Alternative Names ATP synthase subunit d, mitochondrial ATPase subunit d

Function Mitochondrial membrane ATP synthase (F(1)F(0) ATP synthase or Complex

V) produces ATP from ADP in the presence of a proton gradient across the membrane which is generated by electron transport complexes of the respiratory chain. F-type ATPases consist of two structural domains, F(1) - containing the extramembraneous catalytic core, and F(0) - containing the membrane proton channel, linked together by a central stalk and a peripheral stalk. During catalysis, ATP synthesis in the catalytic domain of F(1) is coupled via a rotary mechanism of the central stalk subunits to proton

translocation. Part of the complex F(0) domain and the peripheric stalk, which acts as a stator to hold the catalytic alpha(3)beta(3) subcomplex and subunit

a/ATP6 static relative to the rotary elements.

Cellular Localization Mitochondrion. Mitochondrion inner membrane.

St John's Laboratory Ltd

F +44 (0)207 681 2580

T +44 (0)208 223 3081

W http://www.stjohnslabs.com/ E info@stjohnslabs.com